

IONSCAN™ 500DT

SIMULTANEOUS EXPLOSIVES AND NARCOTICS TRACE DETECTOR



Feature Highlights

- **Simultaneous detection of explosives and narcotics**
- **Over 40 substances detected and identified in 5-8 seconds**
- **Large touch-screen color display**
- **Built-in thermal printer**
- **Internal data storage**

Meets ECAC/EU standard for cargo and passenger screening

The IONSCAN 500DT fills the growing need of security professionals to have the ability to detect a wide range of substances and to be able to adapt as threats and their needs change.

By incorporating two IMS detectors in a single unit, the IONSCAN 500DT is capable of detecting and identifying explosives and narcotics during a single analysis, giving you the ability to detect a broader range of substances while maintaining the high sensitivity and specificity you expect from Smiths Detection.

This dual detector system is ideal for any security application. Aviation security personnel and the military will benefit from the capability to detect a much broader range of explosive substances. Customs, coast guard and border crossing agencies that require the ability to detect both narcotics and explosives from a single sample will be able to execute their duties utilizing one

instrument. Facility security professionals will be aided in better protecting access points to commercial buildings, courthouses and correctional facilities, with the ability to detect multiple threats with one instrument.

Evolution of the IONSCAN

Building on the success and trusted technology of the IONSCAN family of trace detectors, the IONSCAN 500DT has an ergonomic design and a large, color touch-screen display, built-in printer and minimum 300GB internal hard drive. Further advancements include full PC functionality, increased memory and networking capabilities.

The IMS detectors in the IONSCAN 500DT are the same as in the IONSCAN 400B, of which over 14,000 are in use by aviation security, government and law enforcement agencies, the military and by private security companies.

Continued overleaf

IONSCAN 500DT

Sample Collection System

An advanced sampling wand eliminates the need for the operator to handle the sample swab after every sample. This increases the quality of the sample collected and analyzed, further enhancing the performance of the IONSCAN 500DT.

Simplified Maintenance

The IONSCAN 500DT includes a regenerative Air Purification System (APS), which reduces maintenance and cost of ownership.

In addition, most operator replaceable parts do not require the use of tools in the field.

Software

The IONSCAN 500DT utilizes specially developed proprietary software. The operator navigates through the software using the large, color touch-screen display. The large size and color-coded interface of the display makes it easy to view results and operate the IONSCAN 500DT.

With the increased internal memory, more analysis results can be stored. These results can be printed with the built-in printer and exported via the USB port. The software also allows advanced operators the capacity to further review and analyze results.

Technical Data

General Specifications

Technology	Dual Ion Mobility Spectrometry (IMS)
Operating modes	Explosives/Narcotics simultaneous, Explosives only, Narcotics only
Explosives detected	RDX, PETN, NG, TNT, HMX, TATP and others
Narcotics detected	Cocaine, Heroin, Amphetamine, Methamphetamine, MDA, THC and others
Sensitivity	Explosives: picogram range Narcotics: sub-nanogram range
Analysis time	5-8 seconds
Warm-up time	Ready within 30 minutes
Input voltage	95-265 VAC, 50-60Hz 600W cold, 300W warm
Power	110/220V
Weight	19kg (43lbs.)
Display	26.42cm (10.4in), color TFT touch-screen display
Dimensions	40 x 31 x 40cm (16 x 12.5 x 16in) with screen down 40 x 57 x 40cm (16 x 22.5 x 16in) with screen up
Printer	Built-in thermal printer
Operating temperature range	0° to 40°C, (32° to 104°F) less than 95% humidity, non-condensing
Options	Transportation case

For product information, sales or service, please go to www.smithsdetection.com/locations

Smiths Detection, 459 Park Avenue, Bushey, Watford, Herts, WD23 2BW UK
Modifications reserved. 95588905 07/23/15 © Smiths Detection
IONSCAN is a trademark of Smiths Detection Group Ltd.

smiths detection